



Our SC20 and SC20Fast were replaced in 2015 with our market leading SC24 and SC26 field cancelling systems.

At the end of 2024, the SC20 and SC20Fast will be considered End of Life (EOL) by Spicer Consulting and support and repairs of these legacy products will be more restricted.

Upgrade Offers

Upgrade your old SC20 / SC20Fast & get 13% Off new SC24 / SC26 Control Units & Sensors



SC20 → SC24:

- 12 months Warranty with any new items
- Improved user interface
- Auto setup with warning & errors
- Technical support
- AC field cancelling from 0.5 Hz to 5 kHz
- DC field cancelling from DC to 5 kHz
- With AC+DC sensor
 - Cancelling factor at 50/60 Hz increases from 75X to 100X
 - Cancelling factor at DC increases from 200X to 400X

SC20Fast → SC26:

- 12 months Warranty with any new items
- Improved user interface
- Auto setup with warning & errors
- Technical support
- Cancelling factor at 9 kHz increases from 12X to 20X
- Cancelling factor at 50/60Hz increases from 25X to 40X
- Cancelling factor at DC (with DC sensor) increases from 50X to 200X
- Fixes ambient for electron beam tools in 300 mm wafer fabs



	Part Number	Upgrade Price*	Price Saving
SC24 Control Unit	SP124001	€ 14 393	€ 2 151
SC26 Control Unit	SP126001	€ 15 552	€ 2 324
SC24/ DC+ AC Sensor	SP224002	€ 5 592	€ 836
SC24/ AC Sensor	SP224001	€ 2 978	€ 445
SC26/ AC Fast Sensor	SP226001	€ 4 277	€ 639

* Further terms & conditions may apply. Upgrades are one for one basis per existing SC20 or SC20Fast; Only applies to the original purchaser / owner of the SC20 or SC20Fast; A photograph of the existing unit (clearly showing the serial number) will be required; Units may be purchased through a distributor - for Germany and BeNeLux it is point electronic GmbH - or directly from Spicer Consulting; Although in many cases the SC24 and SC26 will be a direct replacement for a SC20/SC20Fast and your existing cables and installation can be utilised, Spicer Consulting cannot be responsible if additional purchases or installation is required to commission a new system. SC26 units will require a new Fast Sensor to maximise performance improvements. Upgrade price represents a 13% off discount on the current end-customer list price. This offer expires on the 31 of December 2024.

Please get in touch for more information.

+49 345 1201190 - www.pointelectronic.de - sales@pointelectronic.de